

Presentation on how to install LAUETOOLSNN

Ravi raj purohit PURUSHOTTAM RAJ PUROHIT¹, Jean-Sebastien MICHA¹

¹French CRG beamline BM32, ESRF, Grenoble

Requirements

- Anaconda distribution
 - <https://www.anaconda.com/products/individual>
- Create a new environment
 - `conda create --name lauenn python=3.8`
 - Lauetoolsnn requires python >3.6
- Activate the created environment and install lauetoolsnn
 - `conda activate lauenn`
 - `pip install lauetoolsnn`
 - Includes all required libraries (in case of errors; try the last command with “--user” suffix)

Home

Environments

Learning

Community



Discover premium data
science content

Documentation

Anaconda Blog



Search Environments



base (root)

bayesian

bm02_beamline

laue

laueGUI

lauetoolsneural



Installed

Channels

Update index...


Search Packages

Newly created environment appears here

Name	T	Description	Version
absl-py		Abseil python common libraries, see https://github.com/abseil/abseil-py.	0.14.0
alabaster		Configurable, python 2+3 compatible sphinx theme.	0.7.12
appdirs		A small python module for determining appropriate platform-specific dirs.	1.4.4
argh		The natural cli.	0.26.2
arrow		Better dates & times for python	0.13.1
astroid		A abstract syntax tree for python with inference support.	2.6.6
astunparse			1.6.3
async-generator			1.10
async_generator		Async generators and context managers for python 3.5+	1.10
atomicwrites		Atomic file writes	1.4.0
attrs		Attrs is the python package that will bring back the joy of writing classes by relieving you from the drudgery of implementing object protocols (aka dunder methods).	21.2.0
autopep8		A tool that automatically formats python code to conform to the pep 8 style guide	1.5.6
babel		Utilities to internationalize and localize python applications	2.9.1
backcall		Specifications for callback functions passed in to an api	0.2.0
bcrypt		Modern password hashing for your software and your servers	3.2.0
binaryornot			0.4.4
black		The uncompromising code formatter	19.10b0
bleach		Easy, whitelist-based html-sanitizing tool	4.0.0
brotlipy		Python bindings to the brotli compression library	0.7.0
build			0.7.0
ca-certificates		Certificates for use with other packages.	2021.9.30
cachetools		Extensible memoizing collections and decorators	4.2.2
certifi		Python package for providing mozilla's ca bundle.	2021.10.8


217 packages available

Conda environment with lauetoolsnn installed



```
C:\Windows\system32\cmd.exe - activate "C:\Users\purushot\Anaconda3\envs\lauetoolsneural" - conda activate lauetoolsneural - lauetoolsnn
(lauetoolsneural) C:\Users\purushot>lauetoolsnn
-- OK! You are using python 3
Missing library libtiff, Please install: pylibtiff if you need open some tiff images
-- warning: module Image or PIL is not installed, but only used for templateimagematching
2021-10-29 15:58:52.207750: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cusolver64_11.dll'; dLError: cusolver64_11.dll not found
2021-10-29 15:58:52.216681: W tensorflow/core/common_runtime/gpu/gpu_device.cc:1835] Cannot dlopen some GPU libraries. Please make sure the missing libraries mentioned above are installed properly if you would like to use GPU. Follow the guide at https://www.tensorflow.org/install/gpu for how to download and setup the required libraries for your platform.
Skipping registering GPU devices...
2021-10-29 15:58:52.220879: I tensorflow/core/platform/cpu_feature_guard.cc:142] This TensorFlow binary is optimized with oneAPI Deep Neural Network Library (oneDNN) to use the following CPU instructions in performance-critical operations:
AVX AVX2
To enable them in other operations, rebuild TensorFlow with the appropriate compiler flags.
Screen: \\.\\DISPLAY1
Size: 3840 x 2160
Available: 3840 x 2080
```

Once installed, launch the GUI by calling **lauetoolsnn** in the terminal



```
python

### config file for LaueNeuralNetwork
## comments

[GLOBAL_DIRECTORY]
prefix =
## directory where all training related data and results will be saved
main_directory = C:\Users\purushot\Desktop\pattern_matching\experimental\GUIv0\latest_version

[MATERIAL]
## same material key as lauetools (see dictlauetools.py for complete key)
## as of now symmetry can be cubic, hexagonal, orthorhombic, tetragonal, trigonal, monoclinic, triclinic

material = In2Bi
symmetry = hexagonal

## if second phase is present, else none
material1 = In_epsilon
symmetry1 = tetragonal

[DETECTOR]
## path to detector calibration file (.det)
detectorfile = C:\Users\purushot\Desktop\In_JSM\calib.det
## Max and Min energy to be used for generating training dataset, as well as for calculating matching rate
emax = 21
emin = 5

[TRAINING]
## classes_with_frequency_to_remove: HKL class with less appearance than specified will be ignored in output
## desired_classes_output : can be all or an integer: to limit the number of output classes
## max_HKL_index : can be auto or integer: Maximum index of HKL to build output classes
## max_nb_grains : Maximum number of grains to simulate per lauepattern
##### Material 0
classes_with_frequency_to_remove = 500
desired_classes_output = all
max_HKL_index = 5
max_nb_grains = 1
##### Material 1
## HKL class with less appearance than specified will be ignored in output
classes_with_frequency_to_remove1 = 500
desired_classes_output1 = all
max_HKL_index1 = 5
max_nb_grains1 = 1

## Max number of simulations per number of grains
## Include single crystal misorientation (1 deg) data in training
## Maximum angular distance to probe (in deg)
## step size in angular distribution to discretize (in deg)
## batch size and epochs for training
max_simulations = 1000
include_small_misorientation = false
angular_distance = 90
step_size = 0.1
batch_size = 50
epochs = 5

[PREDICTION]
# model_weight_file: if none, it will select by default the latest H5 weight file, else provide a specific model
# softmax_threshold_global: thresholding to limit the predicted spots search zone
# mr_threshold_global: thresholding to ignore all matrices less than the MR threshold
# cap_matchrate: any UB matrix providing MR less than this will be ignored
# coeff: should be same as cap_matchrate or no? (this is for try previous UB matrix)
# coeff_overlap: coefficient to limit the overlapping between spots; if more than this, new solution will be computed
# mode_spotCycle: How to cycle through predicted spots (slow or fast or multiorimat) ##slow is more reliable but slow as the name suggests
UB_matrix_to_detect = 1
image_grid_x = 51
image_grid_y = 51

matrix_tolerance = 0.9
matrix_tolerance1 = 0.9
```

Example_config: provides a template to create a config file for automated analysis

Setting in config file

- `### config file for LaueNeuralNetwork`
- `[GLOBAL_DIRECTORY]`
- `prefix =`
- `## directory where all training related data and results will be saved`
- `main_directory = C:\Users\purushot\Desktop\pattern_matching\experimental\GUIv0\latest_version`
-

Prefix: specify a string to distinguish the model files or keep blank

- `[MATERIAL]`
- `## same material key as lauetools (see dictlauetools.py for complete key)`
- `## as of now symmetry can be cubic, hexagonal, orthorhombic, tetragonal, trigonal, monoclinic, triclinic`

- `material = In2Bi`
- `symmetry = hexagonal`
-

If 2 materials present specify them here; in case of only 1 material, keep none for the other

- `## if second phase is present, else none`
- `material1 = In_epsilon`
- `symmetry1 = tetragonal`
-

- `[DETECTOR]`
- `## path to detector calibration file (.det)`
- `detectorfile = C:\Users\purushot\Desktop\In_JSM\calib.det`
- `## Max and Min energy to be used for generating training dataset, as well as for calculating matching rate`
- `emax = 21`
- `emin = 5`
-

Detector settings

Setting in config file

- [TRAINING]
- `## classes_with_frequency_to_remove`: HKL class with less appearance than specified will be ignored in output
- `## desired_classes_output`: can be all or an integer: to limit the number of output classes
- `## max_HKL_index`: can be auto or integer: Maximum index of HKL to build output classes
- `## max_nb_grains`: Maximum number of grains to simulate per lauepattern
- ##### Material 0
- `classes_with_frequency_to_remove = 500`
- `desired_classes_output = all`
- `max_HKL_index = 5`
- `max_nb_grains = 1`
- ##### Material 1
- `## HKL class with less appearance than specified will be ignored in output`
- `classes_with_frequency_to_remove1 = 500`
- `desired_classes_output1 = all`
- `max_HKL_index1 = 5`
- `max_nb_grains1 = 1`
-
- `## Max number of simulations per number of grains`
- `## Include single crystal misorientation (1 deg) data in training`
- `## Maximum angular distance to probe (in deg)`
- `## step size in angular distribution to discretize (in deg)`
- `## batch size and epochs for training`
- `max_simulations = 1000`
- `include_small_misorientation = false`
- `angular_distance = 90`
- `step_size = 0.1`
- `batch_size = 50`
- `epochs = 5`
-

Settings related to Training of Neural network

If unsure, just modify `max_nb_grains` and `max_nb_grains1` ; i.e. maximum grains expected per material

All the other settings can be kept same

Setting in config file

- [PREDICTION]
- # model_weight_file: if none, it will select by default the latest H5 weight file, else provide a specific model
- # softmax_threshold_global: thresholding to limit the predicted spots search zone
- # mr_threshold_global: thresholding to ignore all matrices less than the MR threshold
- # cap_matchrate: any UB matrix providing MR less than this will be ignored
- # coeff: should be same as cap_matchrate or no? (this is for try previous UB matrix)
- # coeff_overlap: coefficient to limit the overlapping between spots; if more than this, new solution will be computed
- # mode_spotCycle: How to cycle through predicted spots (slow or fast or multiorimat) ##slow is more reliable but slow as the name suggests

Settings related to Prediction

- UB_matrix_to_detect = 1

- image_grid_x = 51

- image_grid_y = 51

- matrix_tolerance = 0.9

- matrix_tolerance1 = 0.9

- material0_limit = 1

- material1_limit = 1

- model_weight_file = none

- softmax_threshold_global = 0.85

- mr_threshold_global = 0.80

- cap_matchrate = 0.01

- coeff = 0.3

- coeff_overlap = 0.3

- mode_spotCycle = slow

- ##true for few crystal and preferred texture case, otherwise time consuming; advised for single phase alone

- use_previous = true

-

if unsure, just modify

Grid_x and grid_y : dimension (steps) of scan area

matrix_tolerance and matrix_tolerance1 ; tolerance angle to search for spots: keep 0.5 if unsure

How many grains to detect per material: material0_limit, material1_limit restricts that

use_previous = true; for trying previously identified rotation matrix

mode_spotCycle: graphmode, beamtime mode are the fastest; subsequently on can use also slow mode for more thorough search.

All the other settings can be kept same

Setting in config file

- [EXPERIMENT]
- experiment_directory = C:\Users\purushot\Desktop\In_JSM\ech875_ROI01
- experiment_file_prefix = ech875_ROI01_
- [PEAKSEARCH]
- intensity_threshold = 90
- boxsize = 15
- fit_peaks_gaussian = 1
- FitPixelDev = 15
- NumberMaxofFits = 3000

Settings related to Experiment and misc

material_phase_always_present = 1; forces the neural network to predict always the first grain from material 1; leads to nice plots

- [STRAINCALCULATION]
- strain_compute = true
- tolerance_strain_refinement = 0.7,0.6,0.5,0.4,0.3,0.2
- tolerance_strain_refinement1 = 0.7,0.6,0.5,0.4,0.3,0.2

MTEX output for plotting PF and IPFs ; specify the proper laue group to get proper colors

- [POSTPROCESS]
- hkls_subsets = [1,1,0],[1,0,0],[1,1,1]
- [DEVELOPMENT]
- # could be 1 or 2 / none in case of single phase
- material_phase_always_present = 1
- write_MTEX_file = true
- material0_lauegroup = 9
- material1_lauegroup = 5

Predefined uniform orientation distribution for all symmetries (used during training)

- It is possible to define particular orientations too (in case of preferred texture) → however this will lead to bias in model

